

IN THE UNITED STATES PATENT AND TRADEMARK OFFICE

In re application of

Takumi TERATANI et al. : Group Art Unit: 1633
Serial No. 10/591,407 : Examiner: Quang NGUYEN
Filed: December 8, 2006 :
For: RAT EMBRYONIC STEM

DECLARATION UNDER 37 CFR 1.132

Commissioner for Patents
P.O. Box 1450
Alexandria, Virginia 22313-1450

Sir:

I, Takahiro Ochiya, declare:

That I am a citizen of Japan, and my post office address is National Cancer Center Research Institute, 5-1-1, Tsukiji, Chuo-ku, Tokyo 104-0045 Japan;
That my education and employment history is as shown in my Curriculum Vitae attached hereto (Attachment 1):
That I am a co-inventor of the above-identified U.S. patent application SN 10/591,407 and directed and supervised the following experiment, which was carried out to demonstrate that rat embryonic stem (ES) cells defined in the instant specification cannot be established by the method described in the Examples of Loring (WO 99/27076) (EXPERIMENT 1) and that alkaline phosphatase-positive rat ES cells can be established by the method of claim 8, even when the culture mediums used contain 2% serum (EXPERIMENT 2), the results of which follow hereunder;

[EXPERIMENT 1]

Reproducibility test of rat ES cell establishment
by the method of Loring (WO 99/27076)

[Materials & Methods]

Culture media

1) Medium for ES cell establishment

DMEM High Glucose medium (MP Bio)
2 mM L-Glutamine (MP Bio)
20% Fetal bovine serum (FBS) (Millipore)
1% Non Essential Amino Acid Stock (MP Bio)
1% Antibiotic antimycotic (GIBCO)
0.1 mM 2-Mercaptoethanol(GIBCO)
2000 U/mL mouse leukemia inhibitory factor (mLIF)
(Millipore)

2) Media for maintaining ES cells

i) Maintenance medium A

DMEM High Glucose medium (MP Bio)
2 mM L-Glutamine (MP Bio)
15% FBS (Millipore)
1% Non Essential Amino Acid Stock (MP Bio)
1% Antibiotic antimycotic (GIBCO)
0.1 mM 2-Mercaptoethanol(GIBCO)

ii) Maintenance medium B

DMEM High Glucose medium (MP Bio)
2 mM L-Glutamine (MP Bio)
15% FBS (Millipore)
1% Non Essential Amino Acid Stock (MP Bio)
1% Antibiotic antimycotic (GIBCO)
0.1 mM 2-Mercaptoethanol(GIBCO)
20 ng/mL basic fibroblast growth factor (bFGF) (AbD Serotec)

3) HAT-containing media for maintaining ES cells

i) HAT(+) maintenance medium A

Maintenance medium A
2% HAT Supplement Solution (MP Bio)

- ii) HAT(+) maintenance medium B
Maintenance medium B
2% HAT Supplement Solution (MP Bio)

Feeder cells

Mitomycin C-treated neomycin-resistant mouse embryonic fibroblasts (MEF; Millipore) were used as feeder cells, which were maintained in a medium having the following composition:

- DMEM (GIBCO)
10% FBS (EQUITECH-BIO, Lot No. SFB30-1502)
1% 1x Antibiotic antimycotic (GIBCO)

Mouse ES cells

HPRT-deficient mouse ES cell line ES-E14TG2a was obtained from American Type Culture Collection (ATCC No. CRL-1821).

[Results]

(1) *Determination of culture conditions of mouse ES cells*

To exclude the possibility of contamination of mouse ES cells co-cultured with rat cells, we determined culture conditions that no mouse ES cells survived after the cultivation.

ES-E14TG2a cells were cultured on feeder cells in either of HAT(+) maintenance medium A and HAT(+) maintenance medium B for 3 days. The results are shown in Figure 1. As shown in Figure 1(A) and (B), mouse ES cells died whereas feeder cells survived, when they were co-cultured in the HAT-containing media for 3 days. On the other hand, when they co-cultured in a medium for maintaining ES cells without HAT, colonies of mouse ES cells were observed as shown Figure 1(C).

Therefore, we determined to culture candidates of rat ES cells obtained by co-culture with the mouse ES cells in HAT(+) maintenance medium A or B for 3 days prior to confirmation of the presence of rat ES cell colonies.




Figure 1. Micrographs of mouse ES cells (ES-E14TG2a) cultured on feeder cells in the presence (A, B) and absence (C) of HAT. The mouse ES cells and feeder cells were co-cultured for 3 days in HAT(+) maintenance medium A (A), HAT(+) maintenance medium B (B) and a medium for maintaining ES cells without HAT (C). No mouse ES cell colony was observed due to the death of mouse ES cells, when culturing in HAT-containing media (A, B).

(2) Establishment of rat ES cells

We examined whether or not rat ES cells can be established by co-culturing inner cell mass (ICM) cells obtained from rat blastocysts with mouse ES cells in the same manner as Example 5 of Loring (WO 99/27076).


Uteri of E4.5 to E5 pregnant Wistar rats were perfused with the above-identified Medium for ES cell establishment to give rat blastocysts. Zona pellucida were removed with Tyrode's solution (Ark Resource), and the blastocysts were transferred onto mitomycin C-treated feeder cells in 6-well plates and cultured with the medium for ES cell establishment containing 20% FBS and mLIF. Three days

after, rat ICM cells derived from the blastocysts were excised using Stem Cell Cutting Tool (Vitrolife), incubated in 1 mM EDTA for 30 minutes and subjected to the culture for inducing rat ES cells.

Dispersed ICM cells were sown onto mitomycin C-treated feeder cells in each well of 6-well plates, and mouse ES cells (ES-E14TG2a; 1.5×10^5 cells/well) were added to the each well. These cells were co-cultured in Maintenance medium A (condition 1) or Maintenance medium B (condition 2).

Three days after, ES cell colonies in the wells were treated with 0.25% trypsin/1 mM EDTA for 15 minutes, the dissociated cells were diluted 1:10, transferred onto mitomycin C-treated feeder cells in 6-well plates. The cells were subcultured every 3-5 days for 4 times in the same manner, and then the resulting ES cell colonies were dissociated in the same manner as used in subculture, the dissociated cells were transferred onto mitomycin C-treated feeder cells in 6-well plates without dilution and cultured in the same Maintenance medium except containing HAT (i.e., condition 1: HAT(+) maintenance medium A; condition 2: HAT(+) maintenance medium B) for 3 days, while exchanging the culture medium everyday, to kill the mouse ES cells.



The overall culture protocols are summarized in Figure 2.



(3) Confirmation of establishment of rat ES cells

After the culture in HAT(+) maintenance medium A or B for 3 days for killing the mouse ES cells in (2) above, the presence or absence of ES cell colonies in the wells were observed by microscopy. The results are shown in Figure 3. Under both condition 1 and condition 2, no ES cell colony was observed in any well, demonstrating that rat ES cells cannot be established.

A



B





Figure 3. Micrographs of co-culture of candidates of rat ES cells and mouse ES cells on feeder cells in HAT-containing media. The cells were cultured in HAT(+) maintenance medium A (A) or HAT(+) maintenance medium B (B) for 3 days. Under both conditions, no ES cell colony was observed.

[Conclusion]

It was confirmed that rat ES cells cannot be established by the method described in Example 5 of Loring (WO

[EXPERIMENT 2]

Establishment of rat ES cells using culture media
containing 2% FBS

[Materials & Methods]

Culture media

1) Medium for ES cell establishment

DMEM/F12 1:1 Mixture medium (GIBCO)
2 mM L-Glutamine (MP Bio)
2% Fetal bovine serum (FBS) (Millipore)
1% Non Essential Amino Acid Stock (MP Bio)
1% Antibiotic antimycotic (GIBCO)
1 mM Sodium pyruvate(MP Bio)
0.1 mM 2-Mercaptoethanol(GIBCO)
1% Nucleoside (Millipore)
10% Knockout Serum Replacement (GIBCO)

2) Media for maintaining ES cells

DMEM/F12 1:1 Mixture medium (GIBCO)
2 mM L-Glutamine (MP Bio)
2% Fetal bovine serum (FBS) (Millipore)
1% Non Essential Amino Acid Stock (MP Bio)
1% Antibiotic antimycotic (GIBCO)
1 mM Sodium pyruvate(MP Bio)
0.1 mM 2-Mercaptoethanol(GIBCO)
1% Nucleoside (Millipore)
10% Knockout Serum Replacement (GIBCO)
1000 U/mL rat leukemia inhibitory factor (rLIF)
(Millipore)

Feeder cells

Mitomycin C-treated neomycin-resistant mouse embryonic fibroblasts (MEF; Millipore) were used as feeder cells, which were maintained in a medium having the following composition:

DMEM (GIBCO)

10% FBS (EQUITECH-BIO, Lot No. SFB30-1502)

1% 1×Antibiotic antimycotic (GIBCO)

[Results]

(1) *Establishment of rat ES cells*

We examined whether or not alkaline phosphatase-positive rat ES cells can be established in the same manner as the method described in Example 1 of the instant specification, even when using culture media containing 2% FBS in place of serum-free media.

Uteri of E4.5 to E5 pregnant Wistar rats were perfused with the above-identified Medium for ES cell establishment to give rat blastocysts. Zona pellucida were removed with Tyrode's solution (Ark Resource), and the blastocysts were transferred onto mitomycin C-treated feeder cells in 6-well plates and cultured with the medium for ES cell establishment containing 2% FBS. Seven days after, the outgrowth of blastocysts was mechanically dissociated using Stem Cell Cutting Tool (Vitrolife). The dissociated cells were transferred onto mitomycin C-treated feeder cells in 6-well plates, cultured in the Medium for maintaining ES cells containing 2% FBS for 5 days. The emerged ES cell colonies were dispersed using Stem Cell Cutting Tool (Vitrolife), the dispersed outgrowth was transferred onto mitomycin C-treated feeder cells in 6-well plates, and cultured in the Medium for maintaining ES cells containing 2% FBS for 4 days. The emerged ES cell colonies were dispersed using Stem Cell Cutting Tool (Vitrolife), the dispersed ES cell colonies were transferred onto mitomycin C-treated feeder cells in 6-well plates, and cultured in the Medium for maintaining ES cells containing 2% FBS to establish rat ES cells. The overall culture protocol is summarized in Figure 4.




Figure 4. Culture protocol of rat ES cell establishment

(2) Alkaline phosphatase assay

The alkaline phosphatase activities of the rat ES cells established in (2) above (passage 3) were assayed by staining alkaline phosphatase-positive cells using Vector Blue Alkaline Phosphatase Substrate (Vector Labs). As a result, the rat ES cells established were shown to be alkaline phosphatase-positive (Figure 5). Figure 6 shows a phase-contrast micrograph of the alkaline phosphatase-positive rat ES cell colony A observed in Figure 5.




Figure 5. Staining assay of alkaline phosphatase-positive rat ES cell colonies.




Figure 6. phase-contrast micrograph of the alkaline phosphatase-positive rat ES cell colony A observed in Figure 5.

[Conclusion]

We confirmed that alkaline phosphatase-positive rat ES cells can be established by the method of claim 8 of the instant application, wherein the serum concentration of the culture media is 2%. This demonstrates that the unexpected result (i.e., successful establishment of rat ES cells) commensurate in scope with the claimed invention.

That I further declare that all statements made herein of my own knowledge are true and that all statements made on information and belief are believed to be true; and further that these statements were made with the knowledge that willful false statements and the like so made are punishable by fine or imprisonment, or both, under Section 1001 of Title 18 of the United States Code and that such willful false statements may jeopardize the validity of the application or any patent issuing thereon.

Signed this 18 day of March, 2011.


Takahiro OCHIYA

Curriculum Vitae

Takahiro OCHIYA

Address: c/o National Cancer Center Research Institute, 5-1-1, Tsukiji, Chuo-ku,
Tokyo 104-0045, Japan
Tel: +81-3-3542-2511 (ext. 4800)
Fax: +81-3-5565-0727
E-mail address: tochiya@ncc.go.jp
Title: Chief, Division of Molecular and Cellular Medicine, Translational Research
Group, National Cancer Center Research Institute

Educational History:

1988 Osaka University, Graduate School of Medicine (Ph.D. awarded)

Professional History:

1988-1992 Research Assistant, Institute for Molecular and Cellular Biology, Osaka
University
1993-1998 Section Head, Genetics Division, National Cancer Center Research Institute
1998-2010 Section Head, Section for Studies on Metastasis, National Cancer Center
Research Institute
1991-1992 The Barnham Institute Medical Research, La Jolla, CA USA
2004- Visiting Professor, Waseda University, Graduate School of Science and
Engineering, Bioscience and Biomedical Engineering
2008- Visiting Professor, Tokyo Institute of Technology School and Graduate
School of Bioscience and Biotechnology
2010- Chief, Division of Molecular and Cellular Medicine, Translational Research
Group, National Cancer Center Research Institute

Society: Japan Society of Gene Therapy
Japanese Cancer Association
Japanese Biochemical Society
The Molecular biology Society of Japan

Publication

1. Adenovirus-mediated transfer of the HST-1(FGF4) gene induces increased levels of platelet count in vivo.
Proceedings of the National Academy of Sciences USA (Sakamoto H, Ochiya T, Sato Y, Tsukamoto M, Konishi H, Saito I, Sugimura T, Terada M.) vol.91 No.26, 12368-12372, 1994.
2. Adenovirus-mediated transfer of HST-1/FGF-4 gene protects mice from lethal irradiation.
Oncogene (Y. Takahama, T. Ochiya, H. Tanooka, H. Yamamoto, H. Sakamoto, Hiroshige Nakano and Masaaki Terada.) vol.18 No.43, 5943-5947, 1999.
3. New delivery system for plasmid DNA in vivo using atelocollagen as a carrier material: the Minipellet.
Nature Medicine (T. Ochiya, Y. Takahama, S. Nagahara, Y. Sumita, A. Hisada, H. Itoh, Y.

- Nagai and M. Terada.) vol.6, 707-710, 1999.
4. PDGF alpha-receptor is unresponsive to PDGF-AA in aortic smooth muscle cells from the NG2 knockout mouse.
J Cell Sci. (K. A. Grako, T. Ochiya, D. Barritt, A. Nishiyama and W. B. Stallcup.) vol.112, 905-915, 1999.
5. Detection of spatial localization of Hst-1/Fgf-4 gene expression in brain and testis from adult mice.
Oncogene (Hanako Yamamoto, Takahiro Ochiya, Yasushi Takahama, Yasuo Ishii, Noriko Osumi, Hiromi Sakamoto and Masaaki Terada.) vol.19, 3805-3810, 2000.
6. Mouse flt-1 promoter directs endothelial-specific expression in the embryoid body model of embryogenesis.
Biochem Biophys Res Commun. (G. Quinn, T. Ochiya, M. Terada, and T. Yoshida.) vol.276 No.3, 1089-1099, 2000.
7. Biomaterials for Gene Delivery: Atelocollagen-mediated Controlled Release of Molecular Medicines.
Current Gene Therapy (Takahiro Ochiya, Shunji Nagahara, Akihiko Sano, Hiroshi Itoh, and Masaaki Terada.) vol.1 No.1, 31-52, 2001.
8. High susceptibility of transgenic rats carrying the human c-Ha-ras proto-oncogene to chemically-induced mammary carcinogenesis.
Mutation Research (Tsuda H, Asamoto M, Ochiya T, Toriyama-Baba H, Naito A, Ota T, Sekiya T, Terada M.) vol.477, 173-182, 2001.
9. Involvement of poly(ADP-ribose) polymerase in trophoblastic cell differentiation during tumorigenesis.
Mutation Research (Masutami M, Nozaki T, Watanabe M, Ochiya T, Hasegawa F, Nakagama H, Suzuki H, Sugimura T.) vol.477 No.1-2, 111-117, 2001.
10. Atelocollagen-based gene transfer in cells allows high-throughput screening of gene functions.
Biochem Biophys Res Commun (Honma K, Ochiya T, Nagahara S, Sano A, Yamamoto H, Hirai K, Aso Y, Terada M.) vol.289 No.5, 1075-1081, 2001.
11. HST-1/FGF-4 gene activation induces spermatogenesis and prevents adriamycin-induced testicular toxicity.
Oncogene (Yamamoto H, Ochiya T, Tamamushi S, Toriyama-Baba H, Takahama Y, Hirai K, Sasaki H, Sakamoto H, Saito I, Iwamoto T, Kakizoe T, Terada M.) vol.21 No.6, 899-908, 2002.
12. Enhanced skin carcinogenesis in Cyclin D1-conditional Transgenic mice.
Cancer Research (Yamamoto H, Ochiya T, Takeshita F, Toriyama-Baba H, Hirai K, Sasaki H, Sasaki H, Sakamoto H, Yoshida T, Saito I, Terada M.) vol.62 No.6, 1641-1647, 2002.
13. Mammalian Rcd 1 is a novel transcriptional cofactor that mediates retinoic acid-induced cell differentiation.
EMBO Journal (Hiro N, Ito T, Yamamoto H, Ochiya T, Jinno S, Okayama H.) vol.21 No.19, 5235-5244, 2002.
14. Differentiation of embryonic stem cells into hepatocytes: biological functions and therapeutic application.
Hepatology (Yamamoto H, Quinn G, Asan A, Yamanokuchi H, Taratani T, Terada M, Ochiya T.) vol.37 No.5, 983-993, 2003.
15. Atelocollagen-mediated delivery of synthetic small interfering RNA for effective gene

silencing in vitro and in vivo.

Nucleic Acids Res. (Minakuchi Y, Takeshita F, Kosaka N, Sasaki H, Yamamoto T, Kouno M, Honma K, Nagahara S, Hanai K, Sano A, Kato T, Terada M, Ochiya T.) vol.32, e109, 2004.

16. HST-1/FGF-4 plays critical role in crypt cell survival and facilitates epithelial cell restitution and proliferation.
Oncogene (Sasaki H, Hirai K, Yamamoto Y, Tanooka H, Sakamoto H, Iwamoto T, Takahashi T, Terada M, Ochiya T.) vol.23, 3681-3688, 2004.
17. Direct hepatic fate specification from mouse embryonic stem cells.
Hepatology (Teratani T, Yamamoto H, Aoyagi K, Sasaki H, Asari A, Quinn G, Sasaki H, Terada M, Ochiya T.) vol.41 No.4, 836-846, 2005.
18. Efficient delivery of small interfering RNA to bone-metastatic tumors by using atelocollagen in vivo.
Proceedings of the National Academy of Sciences USA (Takeshita F, Minakuchi Y, Nagahara S, Honma K, Sasaki H, Hirai K, Teratani T, Namatame N, Yamamoto Y, Hanai K, Kato T, Sano A, Ochiya T.) vol.102 No.34, 12177-12182, 2005.
19. Recapitulation of in vivo gene expression during hepatic differentiation from murine embryonic stem cells.
Hepatology (Yamamoto Y, Teratani T, Yamamoto H, Quinn G, Murata S, Ikeda R, Kinoshita K, Matsubara K, Kato T, Ochiya T.) vol.42 No.3, 558-567, 2005.
20. MOZ is essential for maintenance of hematopoietic stem cells.
Genes and Development (Katsumoto T, Aikawa Y, Iwama A, Ueda S, Ichikawa H, Ochiya T, Kitabayashi I.) vol.20 No.10, 1321-1330, 2006.
21. FGF-4 regulates neural progenitor cell proliferation and neuronal differentiation.
FASEB Journal (Kosaka N, Kodama M, Sasaki H, Yamamoto Y, Takeshita F, Takahama Y, Sakamoto H, Kato T, Terada M, Ochiya T.) vol.20 No.9, 1484-1485, 2006.
22. Hedgehog signal activation in gastric pit cell and in diffuse-type gastric cancer.
Gastroenterology (Fukaya M, Isobata N, Ohta H, Aoyagi K, Ochiya T, Saeki N, Yanagihara K, Nakanishi Y, Taniguchi H, Sakamoto H, Shimoda T, Nimura Y, Yoshida T, Sasaki H.) vol.131 No.1, 14-29, 2006.
23. A photon counting technique for quantitatively evaluating progression of peritoneal tumor dissemination.
Cancer Research (Yanagihara K, Takigahira M, Takeshita F, Komatsu T, Nishio K, Hasegawa F, Ochiya T.) vol.66 No.15, 7532-7539, 2006.
24. Adipose tissue-derived mesenchymal stem cells as a source of human hepatocytes.
Hepatology (Banas A, Teratani T, Yamamoto Y, Tokuhara M, Takeshita F, Quinn G, Okochi H, Ochiya T.) vol.46 No.1, 219-228, 2007.
25. Stem cell plasticity: learning from hepatogenic differentiation strategies.
Developmental Dynamics (Banas A, Yamamoto Y, Teratani T, Ochiya T.) vol.236 No.12, 3228-3241, 2007.
26. Atelocollagen for protein and gene delivery.(Review)
Adv Drug Deliv Rev. (Sano A, Maeda M, Nagahara S, Ochiya T, Honma K, Itoh H, Miyata T, Fujioka K.) vol.55 No.12, 1651-1677, 2003.
27. The role of atelocollagen-based cell transfection array in high-throughput screening of gene functions and in drug discovery.(Review)
Curr Drug Discov Technol. (Honma K, Miyata T, Ochiya T.) vol.1 No.4, 287-294, 2004.

28. Therapeutic potential of RNA interference against cancer.(Review)
Cancer Science (Takeshita F, Ochiya T.) vol.97 No.8, 689-696, 2006.
29. "Stem cells into liver"--basic research and potential clinical applications.(Review)
Adv Exp Med Biol. (Banas A, Quinn G, Yamamoto Y, Teratani T, Ochiya T.) vol.585, 3-17, 2006.
30. Application of atelocollagen-mediated siRNA delivery for RNAi therapies.(Review)
Yakugaku Zasshi (Honma K, Takeshita F, Ochiya T.) vol.127 No.5, 807-812, 2007.
31. Atelocollagen-mediated drug discovery technology.
Expert Opin Drug Discov. (Ochiya T, Honma K, Takeshita F, Nagahara S.) vol.2, 159-167, 2007.
32. Establishment of rat embryonic stem cells and making of chimera rats.
PLoS ONE (Ueda S, Kawamata M, Teratani T, Shimizu T, Tamai Y, Ogawa H, Hayashi K, Tsuda H, Ochiya T.) vol.3, issue 7, 2008.
33. A comparative analysis of the transcriptome and signal pathways in hepatic differentiation of human adipose mesenchymal stem cells.
FEBS J (Yamamoto Y, Banas A, Murata S, Ishikawa M, Lim CR, Teratani T, Hatada I, Matsubara K, Kato T, Ochiya T.) vol.275, 1260-1273, 2008.
34. Tissue array substratum composed of histological sections: a new platform for orienting differentiation of embryonic stem cells towards hepatic lineage.
Tissue Eng Part A (Takeuchi T, Ochiya T, Takezawa T.) vol.14, 267-274, 2008.
35. Identification of erythropoietin-induced microRNAs in haematopoietic cells during erythroid differentiation.
Br J Haematol (Kosaka N, Sugiura K, Yamamoto Y, Yoshioka Y, Miyazaki H, Komatsu N, Ochiya T, Kato T.) vol.142, 293-300, 2008.
36. Presentation of functional foreign peptides on the surface of SV40 virus-like particles.
J Biotechnol (Takahashi R, Kanesashi S, Inoue T, Enomoto T, Kawano M, Tsukamoto H, Takeshita F, Imai T, Ochiya T, Kataoka K, Yamaguchi Y, Handa H.) vol.135, 385-392, 2008.
37. RPN2 gene confers docetaxel resistance in breast cancer.
Nat Med (Honma K, Iwao-Koizumi K, Takeshita F, Yamamoto Y, Yoshida T, Nishio K, Nagahara S, Kato K, Ochiya T.) vol.14 No.9, 939-948, 2008.
38. An siRNA against JC virus (JCV) agnoprotein inhibits JCV infection in JCV-producing cells inoculated in nude mice.
Neuropathology.(Matoba T, Orba Y, Suzuki T, Makino Y, Shichinohe H, Kuroda S, Ochiya T, Itoh H, Tanaka S, Nagashima K, Sawa H.) vol.28 No.3, 286-94, 2008
39. Susceptibility of Snark-deficient mice to azoxymethane-induced colorectal tumorigenesis and the formation of aberrant crypt foci.
Cancer Sci.(Tsuchihara K, Ogura T, Fujioka R, Fujii S, Kuga W, Saito M, Ochiya T, Ochiai A, Esumi H.) vol.99 No.4, 677-82, 2008.
40. Down regulation of BRCA2 causes radio-sensitization of human tumor cells in vitro and in vivo.
Cancer Sci.(Yu D, Sekine E, Fujimori A, Ochiya T, Okayasu R.) vol.99 No.4, 810-5, 2008.
41. RNAi-based drug discovery and its application to therapeutics.
IDrugs.(Hokaiwado N, Takeshita F, Banas A, Ochiya T.) vol.11 No.4, 274-8, 2008.
42. Glutathione S-transferase Pi mediates proliferation of androgen-independent prostate

- cancer cells.
Carcinogenesis.(Hokaiwado N, Takeshita F, Naiki-Ito A, Asamoto M, Ochiya T, Shirai T.) vol.29 No.6, 1134-8, 2008.
43. Identification of erythropoietin-induced microRNAs in haematopoietic cells during erythroid differentiation.
Br J Haematol.(Kosaka N, Sugiura K, Yamamoto Y, Yoshioka Y, Miyazaki H, Komatsu N, Ochiya T, Kato T.) vol.142 No.2, 293-300, 2008.
 44. IFATS collection: in vivo therapeutic potential of human adipose tissue mesenchymal stem cells after transplantation into mice with liver injury.
Stem Cells.(Banas A, Teratani T, Yamamoto Y, Tokuhara M, Takeshita F, Osaki M, Kawamata M, Kato T, Okochi H, Ochiya T.) vol.26 No.10, 2705-12, 2008.
 45. Deficiency of Myo18B in mice results in embryonic lethality with cardiac myofibrillar aberrations.
Genes Cells.(Ajima R, Akazawa H, Kodama M, Takeshita F, Otsuka A, Kohno T, Komuro I, Ochiya T, Yokota J.) vol.13 No.10, 987-99, 2008.
 46. Antitumor effect of SN-38-releasing polymeric micelles, NK012, on spontaneous peritoneal metastases from orthotopic gastric cancer in mice compared with irinotecan.
Cancer Res.(Nakajima TE, Yamagihara K, Takigahira M, Yasunaga M, Kato K, Hamaguchi T, Yamada Y, Shimada Y, Mihara K, Ochiya T, Matsumura Y.) vol.68 No.22, 9318-22, 2008.
 47. MicroRNAs as biomarkers and therapeutic drugs in human cancer.
Biomarkers.(Osaki M, Takeshita F, Ochiya T.) vol.13 No.7, 658-70, 2008.
 48. Rapid hepatic fate specification of adipose-derived stem cells and their therapeutic potential for liver failure.
J Gastroenterol Hepatol.(Banas A, Teratani T, Yamamoto Y, Tokuhara M, Takeshita F, Osaki M, Kato T, Okochi H, Ochiya T.) vol.24 No.1, 70-7, 2009.
 49. Pleiotropic function of FGF-4: its role in development and stem cells.
Dev Dyn.(Kosaka N, Sakamoto H, Terada M, Ochiya T.) vol.238 No.2, 265-76.
 50. Screening of potential molecular targets for colorectal cancer therapy.
Int J Gen Med.(Honma K, Takemasa I, Matoba R, Yamamoto Y, Takeshita F, Mori M, Monden M, Matsubara K, Ochiya T.) vol.2, 243-57, 2009.
 51. Dicer is required for maintaining adult pancreas.
PLoS One.(Morita S, Hara A, Kojima I, Horii T, Kimura M, Kitamura T, Ochiya T, Nakanishi K, Matoba R, Matsubara K, Hatada I.) vol.4 issue 1, e4212, 2009.
 52. Local and systemic delivery of siRNAs for oligonucleotide therapy.Methods
Mol Biol.(Takeshita F, Hokaiwado N, Honma K, Banas A, Ochiya T.) vol.487:83-92, 2009.
 53. Vandetanib (ZD6474), an inhibitor of VEGFR and EGFR signalling, as a novel molecular-targeted therapy against cholangiocarcinoma.
Br J Cancer.(Yoshikawa D, Ojima H, Kokubu A, Ochiya T, Kasai S, Hirohashi S, Shibata T.) vol.100 No.8, 1257-66, 2009.
 54. Immunological behavior of enhanced green fluorescent protein (EGFP) as a minor histocompatibility antigen with a special reference to skin isograft and specific regulation of local graft-versus-host reaction (GvHR).
Immunol Lett.(Pan XC, Deng YB, Sugawara Y, Makuuchi M, Okabe M, Ochiya T, Sugiura W, Kitazawa Y, Fuji N, Li XK, Miyamoto M, Kimura H.) vol.123 No.2, 103-13,

- 2009.
55. A mouse model of inducible liver injury caused by tet-on regulated urokinase for studies of hepatocyte transplantation.
Am J Pathol.(Song X, Guo Y, Duo S, Che J, Wu C, Ochiya T, Ding M, Deng H.) vol.175 No.5, 1975-83, 2009.
56. MicroRNA-500 as a potential diagnostic marker for hepatocellular carcinoma.
Biomarkers.(Yamamoto Y, Kosaka N, Tanaka M, Koizumi F, Kanai Y, Mizutani T, Murakami Y, Kuroda M, Miyajima A, Kato T, Ochiya T.) vol.14 No.7, 529-38, 2009.
57. Poly(ADP-ribose) Glycohydrolase deficiency sensitizes mouse ES cells to DNA damaging agents.
Curr Cancer Drug Targets. (Fujihara H, Ogino H, Maeda D, Shirai H, Nozaki T, Kamada N, Jishage K, Tanuma S, Takato T, Ochiya T, Sugimura T, Masutani M.) vol.9 No.8, 953-62, 2009.
58. Mutant mouse p53 transgene elevates the chemical induction of tumors that respond to gene silencing with siRNA.
Cancer Gene Ther.(Tanooka H, Tatsumi K, Tsuji H, Noda Y, Katsube T, Ishii H, Ootsuyama A, Takeshita F, Ochiya T.) vol.17 No.1, 1-10, 2010.
59. Systemic delivery of synthetic microRNA-16 inhibits the growth of metastatic prostate tumors via downregulation of multiple cell-cycle genes.
Mol Ther.(Takeshita F, Patrawala L, Osaki M, Takahashi RU, Yamamoto Y, Kosaka N, Kawamata M, Kelnar K, Bader AG, Brown D, Ochiya T.) vol.18 No.1, 181-7, 2010.
60. Commitment of stem cells into functional hepatocytes.
Differentiation.(Ochiya T, Yamamoto Y, Banas A.) vol.79 No.2, 65-73, 2010.
61. Stem cells for hepatic regeneration: the role of adipose tissue derived mesenchymal stem cells.
Curr Stem Cell Res Ther.(Ishikawa T, Banas A, Hagiwara K, Iwaguro H, Ochiya T.) vol.5 No.2, 182-9, 2010.
62. microRNA as a new immune-regulatory agent in breast milk.
Silence.(Kosaka N, Izumi H, Sekine K, Ochiya T.) vol.1 No.1, 7, 2010.
63. Secretory mechanisms and intercellular transfer of microRNAs in living cells.
J Biol Chem.(Kosaka N, Iguchi H, Yoshioka Y, Takeshita F, Matsuki Y, Ochiya T.) vol.285 No.23, 17442-52, 2010.
64. Combined functional genome survey of therapeutic targets for hepatocellular carcinoma.
Clin Cancer Res.(Satow R, Shitashige M, Kanai Y, Takeshita F, Ojima H, Jigami T, Honda K, Kosuge T, Ochiya T, Hirohashi S, Yamada T.) vol.16 No.9, 2518-28, 2010.
65. Transplantation of a fetal liver cell-loaded hyaluronic acid sponge onto the mesentery recovers a Wilson's disease model rat.
J Biochem. (Katsuda T, Teratani T, Ochiya T, Sakai Y.) vol.148 No.3, 281-8, 2010.
66. Circulating microRNA in body fluid: a new potential biomarker for cancer diagnosis and prognosis.
Cancer Sci. (Kosaka N, Iguchi H, Ochiya T.) vol.101 No.10, 2087-92, 2010.
67. Combined vascular endothelial growth factor-A and fibroblast growth factor 4 gene transfer improves wound healing in diabetic mice.
Genet Vaccines Ther. (Jazwa A, Kucharzewska P, Leja J, Zagorska A, Sierpniowska A, Stepniowski J, Kozakowska M, Taha H, Ochiya T, Derlacz R, Vahakangas E, Yla-Herttuala S, Jozkowicz A, Dulak J.) vol.8, 6, 2010.

68. Enhanced effects of secreted soluble factor preserve better pluripotent state of embryonic stem cell culture in a membrane-based compartmentalized micro-bioreactor. Biomed Microdevices. (Chowdhury MM, Katsuda T, Montagne K, Kimura H, Kojima N, Akutsu H, Ochiya T, Fujii T, Sakai Y.) vol.12 No.6, 1097-105, 2010.
60. Versatile applications of microRNA in anti-cancer drug discovery: from therapeutics to biomarkers. Curr Drug Discov Technol. (Iguchi H, Kosaka N, Ochiya T.) vol.7 No.2, 95-105, 2010.
70. Delivery of small interfering RNA with a synthetic collagen poly(Pro-Hyp-Gly) for gene silencing in vitro and in vivo. Dev Growth Differ. (Adachi T, Kawakami E, Ishimaru N, Ochiya T, Hayashi Y, Ohuchi H, Tanihara M, Tanaka E, Noji S.) vol.52 No.8, 693-9, 2010.
71. Secretory microRNAs as a versatile communication tool. Commun Integr Biol. (Iguchi H, Kosaka N, Ochiya T.) vol.3 No.5, 478-81, 2010.
72. Generation of genetically modified rats from embryonic stem cells. Proc Natl Acad Sci USA, (Kawamata M, Ochiya T.) vol.107, 14223-14228, 2010
73. Gene-manipulated Embryonic Stem Cells for Rat Transgenesis. Cellular and Molecular Life Sciences (Kawamata M and Ochiya T.) in press.